



HAL
open science

Quantification by FTIR spectroscopy of RGII, arabinogalactan-proteins and mannoproteins

Jean Claude J. C. Boulet, Pascale Williams, Thierry Doco, Jean-Michel Roger

► **To cite this version:**

Jean Claude J. C. Boulet, Pascale Williams, Thierry Doco, Jean-Michel Roger. Quantification by FTIR spectroscopy of RGII, arabinogalactan-proteins and mannoproteins. Second International Symposium on Macromolecules and Secondary Metabolites of Grapevine and Wines, Jun 2008, Montpellier, France. 1 p. hal-02756248

HAL Id: hal-02756248

<https://hal.inrae.fr/hal-02756248>

Submitted on 3 Jun 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Quantification by FTIR spectroscopy of RG-II, arabinogalactan-proteins and mannoproteins from wines

Boulet J.C.¹, Williams P.¹, Roger J.M.² and Doco T.¹

¹INRA-UMR1083 Sciences pour l'œnologie, F-34060 Montpellier, France

²CEMAGREF-UMR Information et Technologie pour les Agro-procédés
F-34196 Montpellier, France.

bouletjc@supagro.inra.fr

Polysaccharides are of wide importance in wines. They are involved in many wine properties, such as tartaric and protein stability, filtrability, mellowness, protective colloids, and health safety. So it is important to be able to quantify them. Unfortunately, actually only laboratory methods are available. The most useful of them [1] is based on the quantification of the monomers issued from the polymers after acidic hydrolysis, derivatisation and gas chromatography. This method needs specific material, analytical knowledge and a long analysis time. The aim of this work is to propose an alternative method based on the MIR spectroscopy. This method intends to be sufficiently accurate, easy-to-use and fast, to be a helpful tool in wine processes control such as lees breeding, or addition of pectinolytic enzymes.

Fifteen pure fractions of the 3 main polysaccharides families in wine: Rhamnogalacturonane-II (RG-II), Arabinogalactan-proteins (AGPs) and mannoproteins (MPs) were previously characterised by FTIR. Despite their high spectra collinearity, the differences are strong enough to separate RG-II, AGPs and MPs on a simple principal component analysis (PCA). Three data sets were used for calibration and validation. The first one was obtained by mixing as a Scheffe plan 4 pure fractions: RG-II dimer, AGP0 or AGP-4, and MP0. Two other data sets were obtained from a great variety of wines after a preparation step, in order to eliminate oligosaccharides, polyphenols and wine acids. The spectra were acquired on the powder by FTIR onto a Germanium ATR crystal. The reference analytical values were obtained using the reference chemical method. Each of the RG-II, AGPs or MPs wine concentrations can be calculated with the total polysaccharides weight, and its percent in the powder is given by both reference and FTIR methods.

Previously MPs quantification by FTIR spectroscopy had been proposed [2] then improved [3]. Now this quantification has been extended to the others polysaccharides: AGPs and RG-II. Calibration models were built to give RG-II, AGPs and MPs as a percent of the total polysaccharides. The RMSEP obtained are around 5 %, allowing these models to be used as a process control. But in a few cases we observed a lack of robustness of the FTIR analysis due to the presence of an unknown compound. These cases are easily detected by the spectra shape, because there is a shoulder near 1116 cm^{-1} . We will present possible solutions for the FTIR analysis to take into account this unknown compound.

References:

1. Doco, T., Quéllec, N. & Moutounet, M. (1999). Polysaccharide patterns during the aging of Carignan noir red wines. *American Journal of Enology and Viticulture*, **50**, 25-32.
2. Coimbra, A.C., Barros, A.S., Coelho, E., Gonçalves, F., Rocha, M.S. & Delgadillo, I (2005). Quantification of polymeric mannose in wine extracts by FT-IR spectroscopy and OSC-PLS1 regression. *Carbohydrate Polymers*, **61**, 434-440.
3. Boulet J.C., Doco T. & Roger J.M. (2007). Improvement of calibration models using two successive orthogonal projection methods. Application to quantification of wine polysaccharides. *Chemometrics and Intelligent Laboratory Systems*, **87**, 295-302.