

Efficient screening of odorants' taste dimensions by gas chromatography olfactometry associated taste (GC-OAT) and olfactoscan

Carmen Barba, Elisabeth Guichard, Eliane Lagrange, Noëlle Béno, Thierry

Thomas-Danguin

▶ To cite this version:

Carmen Barba, Elisabeth Guichard, Eliane Lagrange, Noëlle Béno, Thierry Thomas-Danguin. Efficient screening of odorants' taste dimensions by gas chromatography olfactometry associated taste (GC-OAT) and olfactoscan. 11. Wartburg symposium on flavor chemistry & biology, Jun 2016, Eisenach, Germany. 1 p., 2016. hal-02795109

HAL Id: hal-02795109 https://hal.inrae.fr/hal-02795109

Submitted on 5 Jun2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Efficient screening of odorants' taste dimensions by gas chromatography olfactometry associated taste (GC-OAT) and olfactoscan

Carmen Barba^{1,2,3}, Elisabeth Guichard^{1,2,3}, Eliane Lagrange^{1,2,3}, Noelle Beno^{1,2,3} and Thierry Thomas-Danguin ¹ INRA, UMR1324 Centre des Sciences du Goût et de l'Alimentation, F-21000 Dijon, France, ² CNRS, UMR6265 Centre des Sciences du Goût et de l'Alimentation, F-21000 Dijon, France, ³ Université de Bourgogne-Franche Comté, UMR Centre des Sciences du Goût et de l'Alimentation, F-21000 Dijon, France

Introduction

- Many odorants are described with taste attributes [1] although these volatile molecules are not able to activate taste receptors.
- Odor-induced taste is the result of a cognitive process, which depends on individuals' experience and associative memory.
- Those odorants congruent with taste may increase sweet [2] or salty taste perception [3] in foods with respectively less sugar or salt added that is an interesting option for healthier food products formulation.

Materials and methods GAS CHROMATOGRAPHY OLFACTOMETRY ASSOCIATED TO TASTE (GC-OAT) SAFE extract 1μL DB-Wax (30m x 0.32mm x 0.5μm) 12 trained judges Detection Frequency (DF) Test 1: Odor descriptors reen Ioral Figure 1: Chromatogram of multi-fruit juice for test 1 Test 2: Odor associated to taste ALTY Figure 2: Chromatogram of multi-fruit juice for test 2 OLFACTOSCAN (OLFACTOMETER COUPLED TO GC-O) SAFE extract 1uL DB-Wax (30m x 0.32mm x 0.5um 12 trained judges Sweet intensity rating Sweet compounds To rate sweetness by gas chromatography orthonasal olfaction olfactometry Not sweet Verv sweet Olfactoscan is a novel coupled dynamic-dilution system (dynamic-dilution olfactometei olfactometer-gas chromatography olfactometry) that allows on-line evaluation of odor mixture Juice interactions [4].

Water bath 20 °C

GC-OAT analysis is an original and efficient approach for the selection of odor-

Four odorants were found to increase odor sweet perception of multi-fruit juice

(p<0.01). The sweetness enhancement in juice by odor is in the range of 20% to 70%.

inducing-taste compounds that may be used to modulate taste perception in food.

Waste

Figure 3: Olfactoscan

Conclusions

Acknowledgements

Objectives

- 1. To develop a methodology to screen, select and identify odorants associated to taste by GAS CHROMATOGRAPHY OLFACTOMETRY ASSOCIATED TO TASTE (GC-OAT)
- To verify the ability of odorants (previously identified) at two different concentrations to 2. increase sweet odor perception in fruit juice by OLFACTOSCAN

Results

ODORANTS ASSOCIATED TO TASTE IN A JUICE EXTRACT




Figure 4: Aromagram sweet association

Among the 67 olfactive areas detected in an extract of a multi-fruit juice by classical GC-O with a panel of EET, 3 to SALTY, 3 to BITTER and 4 to SOUR by GC-OAT. 12 trained judges, 13 were associated to

Table 1: Identification of odor associated to sweetness compounds

SWEET ENHANCEMENT IN ODOR JUICE

SWEET COMPOUNDS					
N°	DF (%) TEST 1	DF (%) TEST 2	CHEMICAL NAME	ODOR ATTRIBUTES (TEST 1)	TASTE ATTRIBUTES (TEST 2)
1	83	67	n.d.	fruity, sweet	sweet (67%)
2	92	75	methyl 2-methyl butanoate	floral, fruity	sweet (58%) sour (17%)
3	100	75	ethyl butanoate	floral, fruity, sweet, cheese	sweet (50%) sour (17%) bitter (8%)
4	100	75	ethyl 2-methyl-butanoate	fruity, sweet	sweet (58%) bitter (8%) sour (8%)
5	67	75	(E)-β-o-cimene	floral, fruity	sweet (75%)
6	92	75	linalool	floral, fruity, sweet	sweet (67%) sour (8%)
7	75	50	β-damascenone	fruity, sweet	sweet (50%)
8	75	42	phenylmethanol	floral, fruity	sweet (42%)
9	67	67	2-phenylethanol	floral	sweet (58%) salty (8%)
10	58	92	(E)-β-ionone	floral, plastic, solvent	sweet (42%) salty (33%) sour (17%)
11	50	50	n.d.	candy, fruity	sweet (42%) bitter (8%)
12	75	42	furaneol	caramel, sweet	sweet (42%)
13	75	58	γ-decalactone	floral, fruity, sweet	sweet (42%) sour (8%) bitter (8%)




Figure 5: At low concentration 3 compounds: methyl 2-methyl -butanoate, ethyl butanoate and ethyl 2-methylbutanoate are able to increase sweet perception.

Figure 6: At high concentration 4 compounds enhance sweetness in odor juice. Ethyl 2-methyl-butanoate increases sweetness at the two concentrations tested.

INRA

Bibliography

1. Stevenson R. J. and Prescott J, The acquisition of taste properties by odors. Learning and motivation, 26, 433-455, 1995

2. Frank, R. A., and Byrman, Taste-smell interactions are tastant and odorant dependent. J. Chemical Sense, 13, 445-455, 1988

3. Svarifuddin, A., Septier, C., Salles, C. and Thomas-Danguin, T. Reducing salt and fat while maintaining taste: An approach on a model food system. Food Quality and Preference, 48, 59-69, 2016

4. Burseg K. and de Jong C., Application of the olfactoscan method to study the ability of saturated aldehydes in masking the odor of methional. Journal of Agricultural and Food Chemistry, 57, 9086-9090, 2009

cnrs

Carmen Barba has received the support of the EU in the framework of the Marie Skłodowska-Curie H2020-MSCA-IF-2014-655545. We would like to thank Eckes Granini France for providing juice samples and Karine Gourrat from ChemoSens Platform for the technical support.

