


Joanna Fourquet<sup>1,2,3</sup>, Maïna Vienne<sup>1,2\*</sup>, Jean Mainguy<sup>1,2\*</sup>, Vincent Darbot<sup>3\*</sup>, Pierre Martin<sup>1,2</sup>, Olivier Bouchez<sup>4</sup>, Adrien Castinel<sup>4</sup>, Sylvie Combes<sup>3</sup>, Carole Iampietro<sup>4</sup>, Christine Gaspin<sup>1,2</sup>, Denis Milan<sup>4</sup>, Cécile Donnadieu<sup>4</sup>, Céline Noïrot<sup>1,2</sup>, Geraldine Pascal<sup>3</sup> and Claire Hoede<sup>1,2\*</sup>

- 1 Université Fédérale de Toulouse, INRAE, BioinfOmics, GenoToul Bioinformatics facility, 31326, Castanet-Tolosan, France  
2 Université Fédérale de Toulouse, INRAE, MIAT, 31326, Castanet-Tolosan, France  
3 GenPhySE, Université de Toulouse, INRAE, ENVT, F-31326, Castanet Tolosan, France  
4 INRAE, GeT-PlaGe, Genotoul – INRAE – 31326 Castanet-Tolosan, France  
\* Present at ECCB 2022

Corresponding author: [claire.hoede@inrae.fr](mailto:claire.hoede@inrae.fr)

## Production of whole metagenome assembly, functional and taxonomic profile




## Workflow features



### Type of NGS data:

whole genome shotgun sequencing (Illumina HiSeq3000 or NovaSeq, paired, 2\*150bp ; PacBio HiFi reads, single-end)




### Workflow:

a scalable and reproducible metagenomic analysis with a nextflow [1] pipeline using Singularity [21] containers




### Fully documented

<https://forgemia.inra.fr/genotoul-bioinfo/metagwgs>



The strategy of aligning reads of every samples (Cross-alignment: All) against each assembly improved the quality of binning. Tests were done on a synthetic mock composed of 142 bacteria and archea genomes with 3 samples of 66.651.100 Illumina paired-end reads (2x150bp).



Execution time in minutes of the original MetaWRAP bin refinement module [25] compare to the improved version implemented in metagWGS, on the synthetic mock. The improved version uses Checkm2 [28] instead of Checkm1 [29] and takes advantage of a custom resume parameter. The modified version gives very similar results.

## Conclusion

The new version of metagWGS (2.3) allows the analysis of Illumina short reads or PacBio HiFi long-reads sequencing data and brings as a major new feature the binning of contigs.

The workflow proposes to use the abundance information contained in nearby samples to improve the binning by implementing cross-alignment per sample set.

We have also improved the performance of the bins refinement step by dividing the execution time by 7.

## References



## Acknowledgements

Projet cofinancé par le Fonds Européen de Développement Régional

ATB\_Biofilm funded by PIREST Anses, 2020/01/142  
Funded by France Génomique (ANR-10-INBS-09-08)  
Funded by Antiselfish Project (Labex Ecofect)  
Funded by ExpoMicoPig project (France Futur élevage)  
Vincent Darbot's work has been supported by RESALAB OUEST  
SeqOcln financed by FEDER funds  
(Programme Opérationnel FEDER-FSE\_Midi-Pyrénées et Garonne 2014-2020)



## Perspectives

Better assembly of the minority species when the sequencing depth is not sufficient: implementation of co-assembly (by giving the possibility to normalize data first).

Improve the performance of the workflow: replacing Prokka [15] with other tools.

Long term perspectives: enable the annotation of antibiotic resistance genes and of the mobileome