

Characterisation of the key aroma compounds in alcohol free beer base by aroma extract dilution analysis

José A. Piornos, Dimitrios P Balagiannis, Elisabeth Koussissi, Eric Brouwer, Jane K Parker

▶ To cite this version:

José A. Piornos, Dimitrios P Balagiannis, Elisabeth Koussissi, Eric Brouwer, Jane K Parker. Characterisation of the key aroma compounds in alcohol free beer base by aroma extract dilution analysis. 15th Weurman Flavour Research Symposium, Sep 2017, Graz, Austria. hal-03839579

HAL Id: hal-03839579 https://hal.inrae.fr/hal-03839579

Submitted on 4 Nov 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Characterisation of the key aroma compounds in alcohol free beer base by aroma extract dilution analysis



HEINEKEN



José A. Piornos¹ | Dimitrios P. Balagiannis¹ | Elisabeth Koussissi² | Eric Brouwer² | Jane K. Parker¹

- ¹ The Flavour Centre, Department of Food and Nutritional Sciences, UK.
- ² Heineken Supply Chain BV, Global Innovation & Research, Zoeterwoude, The Netherlands.

Alcohol Free Beer Cold Contact Process







Mashing Water is added and starch is broken down into fermentable sugars by amylases



Wort boiling At this step hops are added and the wort is boiled.



Centrifuging and cooling Spent grains and hops are separated from the liquid and this is cooled down.



Cold Contact Fermentation Wort is fermented at low temperature (\sim 0°C) to prevent ethanol formation.



Post-treatments Maturation, flavour addition

and bottling, among others.

Introduction

Alcohol-free beers brewed by cold contact process usually exhibit a flavour reminiscent of wort and lack fruitiness [1]. Although the addition of flavourings after fermentation is a common practice, the typical "worty" aroma cannot easily be masked. The aim of this research is to identify the compounds contributing to the characteristic aroma of these beers by means of the Sensomics approach.

Materials and methods

Base beer sample

Brewed by Cold Contact

Alcohol free product (<0.05%) NO external flavours

Process


added to "Base" beer Pasteurised and bottled




Two fibres-HS-SPME


Volatile compounds were adsorbed/absorbed onto the SPME fibre. Simultaneous extraction of two fibres:







Isolation of volatile compound fractions by SAFE^[2]




Identification of odour-active compounds


GC-MS and GC-Olfactometry were performed on SAFE and SPME extracts.

Criteria for identification:

- LRI on polar (FFAP) and non-polar (DB-5) columns
- Odour description
- Bibliography
- Agreement with analyses of authentic compounds




Aroma extract dilution analysis (AEDA)



FD factor for a specific odour region is defined as the highest dilution at which the odour can still be perceived at the sniffing port of the GC-O.

Results and discussion





- 22 odour-active compounds had $FD \ge 16$; with methional, butanedione, 4-vinylphenol, and 2-methoxy-4-methyphenol amongst the highest.
- 4 additional compounds were found in **SPME-GC-O** experiments.
- No fruity ester, typically found in lager and other beers [2], was among these compounds, explaining the lack of fruitiness.
- Strecker aldehydes (methylpropanal, 2-methylbutanal, 3methylbutanal, phenylacetaldehyde, methional) have been reported as contributors to the malty, worty aroma of alcohol free beers [1].
- Identification of 2-acetyl-1-pyrroline is not conclusive.
- Furthermore, 2 compounds remain unidentified.

Conclusions

- The lack of fruity flavour was associated with the absence of odour-active fruity esters.
- Methional was a potential contributor to the "worty" aroma of the alcohol free base beer studied.
- Quantification of these compounds is of great importance to understand their role in the aroma of the product.

References

- 1. Perpète, P. & Collin, S. J. Agric. Food Chem., 1999, 47 2374-2378.
- 2. Langos, D., Granvogl, M., & Schieberle, P. J. Agric. Food Chem., **2013**, 61, 11303-11311.

Acknowledgements

- This project is sponsored by Heineken Supply Chain, Global Innovation & Research.
- Icons retrieved from www.flaticon.com, authored by Freepik (snowflake, peach, flower, honeycomb), Smashicons (nose) and Madebyoliver (cheese, potatoes).

Contact information

- Department of Food and Nutritional Sciences, University of Reading, Whiteknights, RG6 6AH
- Email: j.k.parker@reading.ac.uk | www.reading.ac.uk/flavourcentre/