

MS2Snoop: an R package for automatic extraction and curation of MS/MS spectra in metabolomics

Kevin Wagner, Lain Pavot, Franck Giacomoni, Guillaume Marti, Yann Guitton, Laurent Debrauwer, Emilien L. Jamin, Jean-François Martin

▶ To cite this version:

Kevin Wagner, Lain Pavot, Franck Giacomoni, Guillaume Marti, Yann Guitton, et al.. MS2Snoop: an R package for automatic extraction and curation of MS/MS spectra in metabolomics. Analytics 2022, Sep 2022, Nantes, France. 2022. hal-04001175

HAL Id: hal-04001175 https://hal.inrae.fr/hal-04001175

Submitted on 22 Feb 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



MS2Snoop: a R package for automatic extraction and curation of MS/MS spectra in metabolomics







Kévin Wagner^{1,2}, Lain Pavot³, Franck Giacomoni³, Guillaume Marti^{2,4}, Yann Guitton⁵, Laurent Debrauwer^{1,2}, <u>Emilien Jamin</u>^{1,2}, Jean-François Martin^{1,2}

¹ Toxalim (Research Center in Food Toxicology), Toulouse University, INRAE, ENVT, INP-Purpan, Toulouse, France
 ² MetaboHUB-MetaToul, National Infrastructure of Metabolomics and Fluxomics, Toulouse, France
 ³ Université Clermont Auvergne, INRAE, UNH, Plateforme d'Exploration du Métabolisme, MetaboHUB Clermont, Clermont-Ferrand, France
 ⁴ Metatoul-AgromiX Platform, LRSV, Université de Toulouse, CNRS, UT3, INP, Toulouse, France
 ⁵ Oniris, INRAE, LABERCA, Nantes, France

Introduction

Omics approaches including metabolomics generate huge datasets requiring powerful processing tools. In the field of non-targeted mass spectrometry metabolomics, HRMS data mining has been the subject of much work, with the development of powerful interactive tools, including Workflow4Metabolomics (W4M) [1-2]. One of the main bottlenecks of current workflows lies in the complexity of the interpretation of results, in particular a lack of reliable MS/MS spectral reference databases (DBs) to be used for confident identification (level 1 and 2 according to the metabolomics standard initiative). Creating reliable and well curated spectral DBs remains a laborious and challenging task, but necessary to ensure the quality of the reference spectra. In this work, we present a tool inspired by expert manual curation allowing the export of MS/MS data from both targeted and data dependent (DDA) fragmentation methods to laboratory spectral DBs, named "MS2Snoop".

MS2snoop Workflow (R): **XCMS** frag4feature MS2snoop **msPurity FilterFrag Peak Fragment ions Filtration Precursors Datasets** MS/MS spectra picking reconstruction quality merging Definition of detection level of fragment ions (intensity, From mzML file(s) containing From same mzML file(s) Association of fragment ions Merging of fragment & S/N ratio) MS & MS/MS Use of purityA function to precursor ions with list of with precursor by merging Generation of .tsv files for Use of xcmsSet & group expected compounds evaluate the purity of « purityA » and « xcmsSet »

objects

Examples of applications

Generation of an « xcmsSet »

functions of XCMS


object (.tsv)

Targeted MS/MS by resonant CID achieved on a LTQ-Orbitrap XL using ESI in the positive mode

precursors [3]

object (.tsv)

Generation of a « purityA »




precursor ion intensity & fragment ion intensities

m/z 126 kept by MS2snoop

Curation by correlation between

Generation of .tsv files for

fragment & precursor ions




Successful application to the particular case of acetochlor with a very weak precursor ion intensity

Correlation of fragment

intensity

intensity with precursor

fragment & precursor ions



Conclusions

MS2Snoop already allows extracting and cleaning MS/MS datasets using commonly used metabolomics data processing steps as a starting point, then combined with inputs of a list of expected compounds and MS/MS raw data deconvoluted by MSPurity. Firstly developed and validated with targeted MS/MS spectra from about 30 compounds, MS2Snoop is compatible with DDA spectra, although limited by the number of MS/MS scans acquired. This package will allow high throughput extraction of MS/MS spectra in order to feed the PeakForest project [4]. Preliminary results obtained with MS2Snoop allowed considering it as a valuable tool for confident annotation of metabolites on chromatographic methods coupled to non-targeted mass spectrometry.

A beta version of this original package is already publicly available on W4M (https://workflow4metabolomics.org/) and recent update allows fragments ions formula annotation with Sirius.

Acknowledgments

HRMS analyses were achieved in the frame of the French national infrastructure MetaboHUB (MetaboHUB-ANR-11-INBS-0010). Special thanks to JF Martin for his long career and the good times!



[1] Giacomoni F. et al. Bioinformatics, (2015) 31, 1493 [2] Guitton Y. et al, Int. J. Biochem. Cell Biol., (2017) 93, 89

[3] Lawson T. N. et al., Anal. Chem., (2017) 89, 1432 [4] Paulhe N. et al., Metabolomics, (2022)18, 40

Poster presented at Analytics2022, Nantes, France