



**HAL**  
open science

# Characterization of the genetic and functional diversity of the NLRome/resistome in melon (*Cucumis melo* L.) using Nanopore Adaptive sampling (NAS)

Javier Belinchon-Moreno

► **To cite this version:**

Javier Belinchon-Moreno. Characterization of the genetic and functional diversity of the NLRome/resistome in melon (*Cucumis melo* L.) using Nanopore Adaptive sampling (NAS). Journées Jeunes Chercheurs 2023, INRAE GDEC, Apr 2023, Clermont - Ferrand, France. hal-04071831

**HAL Id: hal-04071831**

**<https://hal.inrae.fr/hal-04071831v1>**

Submitted on 17 Apr 2023

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial - NoDerivatives 4.0 International License

➤ Exploring the genetic and functional diversity of the NLRome/resistome in melon (*Cucumis melo* L.) using Nanopore Adaptive sampling (NAS)

Javier BELINCHON-MORENO

Journées Jeunes Chercheurs, Clermont-Ferrand, 3rd April 2023




# ➤ Oxford Nanopore Sequencing Technology (ONT)

Nanopore embedded in an electrically resistant polymer membrane

Read from *cis* to *trans*

Constant voltage applied → Read change of voltage




Flongle	MinION	MinION <sub>Mk1C</sub>	GridION	PromethION <sup>24</sup>	PromethION <sup>48</sup>
					
Theoretical Max: 2.8 Gb	Theoretical Max: 50 Gb	Theoretical Max: 50 Gb	Theoretical Max: 50 Gb x 5 = 250 Gb	Theoretical Max: 290 Gb x 24 = 7 Tb	Theoretical Max: 290 Gb x 48 = 14 Tb

(Nanopore London Callings, 2022)

## ➤ What is Nanopore adaptive sampling?


- Accept/reject molecules based on **small initial part** of sequence (~500 bp)
- **Enrich/deplete** specific regions of interest.
- Increase on-target data, reduce time-to-answer.
- **Real time** target selection: **easy** library preparation.



(Nanopore London Callings, 2022)

# ➤ Our objective: NAS to decipher the genetic and functional diversity of the NLRome in melon

- 1 Reference genome not enough → Construction of **pan-NLRome**  
 High number of varieties (~150 varieties) → Need for **selective sequencing**



2

Understand its role in the expression of immunity against a wide range of pests and pathogens (GWAS)

(Adapted from Barragan & Weigel, 2021)

# ➤ NLR resistance genes: a conserved structure

**Nucleotide binding site (NB) -leucine rich repeat (LRR)**



**Intracellular** immune receptor proteins

**Major family** of plant-resistance (R) genes

**Broad range of resistance**



*Adapted from (Gottin et al., 2021)*



*F. oxysporum*



*A. gossypii*



Oidium



*P. infestans*


Usually grouped into **clusters**

**High level of presence/absence polymorphisms (PAV)**

Low frequency in the Cucurbitaceae family

# ➤ Our target regions: 15 NLR-gene clusters


Reference genome → Variety **Anso77** (draft genome produced at INRAE-GAFL/EPGV)




Reference genome (.fasta)  
+  
Target regions (.bed)

Clusters:  
NLR genes + 20kb-extra flanking  
Variable in size  
Variable in number of genes

Genome length ~380Mb  
Target regions length ~6.16Mb



# ➤ Our first experience: Variety Anso-77 (same as provided reference)





NAS

Standard sequencing





**Not** masking repetitive elements

**Masking** repetitive elements >200pb






# ➤ Visual enrichment of our target regions



— ROI  $\pm 20\text{kb}$   
— ROI delimiter

## ➤ Good enrichment of our target regions




$$Enrichment_{region} = \frac{\frac{region\_cov\_NAS}{chrom\_cov\_NAS}}{\frac{region\_cov\_standard}{chrom\_cov\_standard}}$$

$$Av. enrichment = \sum_{region=1}^N \frac{enrichment_{region}(X) * region\_size(bp)}{total\_size\_regions(bp)}$$

- Variable between clusters
- Very high in average

$$\frac{region\_cov\_NAS}{region\_cov\_standard} \rightarrow \sim 3-4$$


# ➤ Genome assemblies: Good assembly of target regions compared to a standard sequencing



## Assembly using WGS data:

- Target regions **fragmented** into multiple contigs
- Target regions **poorly covered** compared to adaptive sampling

# ➤ Genome assemblies: Good assembly of target regions




ANSO-77 – inodorus (ssp. melo)

Target regions fully assembled

All predicted NLR genes (84) found in the assembly

number_contigs	15
total_length	7.023 Mb
largest_contig	1.443 Mb
GC(%)	33.23
N50	1.107 Mb
N75	1.026 Mb
L50	3
L75	5

# ➤ New experiences: Sequencing of more genetically distant varieties from Anso77



Single variety on promethION



2 barcodes on promethION




4 barcodes on promethION



8 barcodes on promethION

# ➤ Preliminary results using different varieties: Doublon

Doublon: predicted 15 gene clusters



Doublon assembly of regions of interest

**DOUBLON – cantalupensis (ssp. melo)**


Target regions fully assembled  
**All predicted NLR genes (76) found in the assembly**

number_contigs	15
total_length	7.131 Mb
largest_contig	1.470 Mb
GC(%)	33.31
N50	1.060 Mb
N75	1.003 Mb
L50	3
L75	5



# ➤ Preliminary results using different varieties: Zhimali

Anso-77: predicted 15 gene clusters



Zhimali assembly of regions of interest

ZHIMALI – chinensis (ssp. agrestis)

Target regions fully assembled

**78 NLR genes found in the assembly**

number_contigs	15
total_length	7.635 Mb
largest_contig	1.476 Mb
GC(%)	33.23
N50	1.086 Mb
N75	1.035 Mb
L50	4
L75	5



# ➤ Perspectives

## Short-term perspectives:

1. **Increase gradually the number of barcodes** on PromethION.
2. Start evaluating the available methods for pangenome construction.

## Middle-term perspectives:

1. **Find the optimal number of barcodes** on PromethION: Selective sequencing of 150 varieties.
2. Evaluate the ability to detect **structural variants** and **copy number variations** between varieties.

## Long-term perspectives:

1. Perform a similar experience using a **diploid heterozygote** species with a **higher number of NLR genes** and **more complex clusters: *Vitis vinifera***.



## ➤ In conclusion

### What Nanopore adaptive sampling offers?

- Nanopore adaptive sampling is a **simple, reliable, efficient** and **cost-saving** approach.
  - ✓ Standard DNA extraction (no need for HMW DNA)
  - ✓ Standard library preparation
  - ✓ No need of probes design/synthesis
  - ✓ No need of short-read polishing with Kit14
  - ✓ Real-time sequencing
  - ✓ Reduced off-target volume of data compared to WGS

### What can we do?

1. Adaptive sampling allows to efficiently retrieve (**mapping** and **assembly**) our **15 ROI** in melon.
2. It will allow to **explore the genetic and functional diversity** of melon **NLR genes** using a **large amount of varieties**.

# Acknowledgements



## INRAE-EPGV

Patricia Faivre-Rampant  
Damien Hinsinger  
Aurelie Berard  
Aurelie Canaguier  
Isabelle Le-Clainche

## INRAE-GAFL

Nathalie Boissot  
Veronique Chovelon  
Jacques Lagnel  
Vincent Rittener-Ruff

## Genoscope

Stefan Engelen  
Corinne Cruaud

